ENGLISH

Before beginning your treatments please review this important information.

1. PRODUCT DESCRIPTION

a. What is it, and how does it work?

SkinPen® Precision System is a minimally invasive, micro-needling device intended for use on facial acne scars. The SkinPen Precision System is comprised of a reusable motor unit designed to be attached to sterile, disposable cartridges that house 14 micro- needles (~250 um diameter) on a reciprocating head which when activated and placed properly against the skin can create hundreds to thousands of “micro” punctures into the skin. The procedure creates many microscopic punctures in the skin’s outer layers. As these punctures heal, the remodeling process works to improve the appearance of facial acne scars. SkinPen Precision’s patent-pending cartridge design inhibits the entrance of fluids into the device which along with the BioSheath and routine cleaning help to prevent potential cross-contamination of the motor unit between uses. It is for prescription use only.

b. Who is a candidate?

The SkinPen Precision is intended to utilize Collagen Induction Therapy to treat Melasma, acne and surgical scars. The SkinPen® Precision System’s clinical study was conducted on adult men and women ranging from light to very dark skin (Fitzpatrick skin types).

2. CONTRAINDICATIONS

a. Are there any reasons why I should not receive a SkinPen® Precision treatment?

The use of the SkinPen® Precision System should not be used on patients who:

- Have active skin cancer in the treatment area(s)
- Have open wounds, sores, or irritated skin in the treatment area(s)
- Have an allergy to stainless steel or anesthetics
- Have a hemorrhagic (bleeding) disorder or hemostatic (bleeding) dysfunction
- Are pregnant or nursing
- Are currently taking drugs with the ingredient isotretinoin (such as Accutane)

NOTE: This product is not intended for transdermal (under the skin) delivery of topical products such as cosmetics, drugs, or biologics.

3. PRECAUTIONS

a. What precautions should my doctor advise me about?

The SkinPen® Precision System has not been evaluated in the following patient populations. Therefore, if you have a history of the following conditions or have taken the following medications, please let your doctor know, as treatment with the SkinPen Precision System may not be appropriate for you: Actinic (solar) keratosis; active acne; collagen vascular diseases or cardiac abnormalities; diabetes; eczema, psoriasis and other chronic conditions in the treatment area or on other areas of the body; immunosuppressive therapy; history of contact dermatitis; raised moles in the treatment area; rosacea; active bacterial or fungal infection, active viral herpes simplex infections (such as cold sores); warts; keloid scars; patients on anticoagulants (also known as ‘blood thinners’); scars and stretchmarks less than one year old; scleroderma; and wound-healing deficiencies. Inform your provider if you have any of the above conditions. Your doctor will determine whether treatment with the SkinPen Precision System is right for you.

Let your doctor know if you are allergic or sensitive to any of the following ingredients which are in the Skinfuse Lift HG hydrogel: purified water, glycerin, carbomer, potassium hydroxide, disodium EDTA, phenoxythanol, caprylyl glycol, sorbic acid. If you are, SkinPen® Precision treatment may not be safe for you.

4. BEFORE TREATMENT INFORMATION

a. What happens in the office before the SkinPen® Precision treatment?

Note that each doctor may have a unique process for assessing and treating patients. The following is an example of what you would likely experience with a typical SkinPen® Precision treatment. Before the SkinPen® Precision treatment, your doctor will ask you questions about your medical history, as well as your treatment goals. Your doctor will discuss whether you are an appropriate candidate for a SkinPen® Precision System treatment and review what to expect during and after a treatment, including common treatment responses and adverse events.

You may be cautioned to avoid sun exposure and stop topical retinoid therapy 24 hours prior to procedure. You should allow at least 24 hours after autoimmune therapies before a treatment. You should wait six months following oral isotretinoin (such as Accutane) use before receiving treatment.

Your doctor will also examine your treatment area, and may take photos. Different options for pain management will be discussed, and if pretreatment numbing is desired, a topical anesthetic agent may be used. The treatment area will be cleaned and then prepared with isopropyl alcohol or other antiseptic before the treatment.

5. TREATMENT DESCRIPTION

a. What happens during the treatment?

A layer of Skinfuse Lift HG is applied to the treatment area to protect the skin against abrasion and friction of the SkinPen® Precision device during treatment. **NOTE:** Let your doctor know if you are allergic or sensitive to any of the following ingredients which are in the Skinfuse Lift HG hydrogel: purified water, glycerin, carbomer, potassium hydroxide, disodium EDTA, phenoxythanol, caprylyl glycol, sorbic acid. If you are, the SkinPen® Precision treatment may not be safe for you. The doctor then selects a depth on the microneedle cartridge to begin treatment on the patient and begins the treatment. As the SkinPen® Precision glides over the skin micro-channels are created on the skins surface. Gauze may be used to pat down the affected area after treatment, and it is suggested to apply a generous layer of Skinfuse Lift HG to prevent the skin from drying out post treatment.

b. Does microneedling with the SkinPen® Precision hurt?

The SkinPen® Precision treatment may cause some minor discomfort during the treatment. Your doctor may recommend that a topical numbing agent be used during treatment to further minimize discomfort.

6. AFTER TREATMENT INFORMATION

a. What should I expect following the treatment?

In the SkinPen® Precision clinical study, the most common treatment responses experienced were:

Dryness, rough skin, tightness, redness, itching, peeling, discomfort, tenderness, and burning. These conditions resolved over time without any further complications.

Although not seen in the clinical study, you may experience reactivation of herpes simplex virus (cold sore), pigment changes that include lighter or darker skin in the area treatment that resolves over time, or no change in their acne scars.

Your doctor will also tell you what to expect following a SkinPen® Precision system treatment. Within the first 72 hours post-treatment you should avoid sweaty exercise and sun exposure. Exposure to these conditions could lead to: itching, burning, stinging, and tingling), scaling/dryness, redness, swelling, and tenderness/discomfort.

b. Will I need more than one treatment to achieve my desired results?

You should discuss treatment goals with your doctor. In the clinical study patients were treated in a series of 3 treatments spaced 4 weeks apart. It is recommended to avoid other facial aesthetic treatments the month following your SkinPen Precision treatment.

c. Do the results last forever?

While individual results may vary, in the clinical study, the results were evaluated 6 months after the third (final) treatment. After this, additional treatments may be needed to maintain your desired result.

7. WHEN TO CALL YOUR DOCTOR

a. When should I call my doctor?

Call your doctor immediately if you experience:

1. Signs of infection such as an increase in redness, warmth, itching, or pus formation. This would typically happen within a day or two after the treatment.
2. Allergic reactions to the test material (s). Rare allergic reactions can consist of severe contact dermatitis (rash) or itching, swelling, or difficulty breathing. This risk is increased for individuals with a history of allergies, and individuals with asthma and/or a history of hives, and itching.

Be sure to call your doctor if you:

1. Any additional treatment responses not discussed above.

8. ADDITIONAL INFORMATION

a. What should I do if I have additional questions?

For further questions and information, please contact your SkinPen Treatment Provider.

b. Lift info:

Lift HG is a hydrogel wound dressing (without drugs and/or biologics) to protect against abrasion and friction during the microneedling procedure. It may be applied additionally the day of the procedure to prevent the skin from drying out post procedure.

DIRECTIONS: AM ☼ PM ☾, Apply a generous amount once or twice daily on irritated, sensitive or compromised skin. Reapply as needed.

WARNING: Do not use in the eyes. Keep out of reach of children.

Non-Cytotoxic, Non-Sterile, 0.51 fl oz/15mL

INGREDIENTS: Purified water (Aqua), Glycerin, Carbomer, Potassium Hydroxide, Disodium EDTA, Phenoxythanol, Caprylyl Glycol, Sorbic Acid.

c. Rescue info:

DIRECTIONS: AM ☼ PM ☾, Apply the cosmetic calming complex liberally as needed to help soothe and calm the skin. For use after 24 hours post-treatment. 0.34 fl oz/10mL

INGREDIENTS: Water (Aqua), Glycerin, Hydroxyethylcellulose, Magnesium Aspartate, Zinc Gluconate, Copper Gluconate, Colloidal Oatmeal, Aloe Barbadensis Leaf Juice, Sodium Hyaluronate, Allantoin, Bisabolol, Hydroxymethoxyphenyl Decanone, Niacinamide, Lysophosphatidic acid, Lecithin, Lysolecithin, Acetyl Carboxymethyl Cocoyl Glycine, Ethyl Macadamiate, Borago Officialis Seed Oil, Tocopheryl Linoleate/Oleate, Xanthan Gum, Maceratum Gum, Caprylyl Glycol, Phenoxyethanol, Tromethamine.

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